Analyzing current vs voltage relation in resistors, potentiometers and lightbulbs

Sebastien Psarianos Sofiya P'yavka

November 1, 2022

1 Introduction

The objective of this experiment was to determine the charge to mass ratio of an electron. An electron moving in a magnetic field will experience a force of:

$$F = e\vec{\nu} \times \vec{B}$$

e is the charge of the electron, v is its velocity and B is the magnetic field. Using Newton's second law of motion, this can be further written in terms of $F = ma = m\frac{v^2}{r} = evB$ where r is the radius of the circular orbit that the electron will follow and m is its mass. The electron will be accelerated through a potential difference in this experiment and the curvature of the path can be described using:

$$\frac{1}{r} = \frac{B}{\sqrt{V}} \sqrt{\frac{e}{2m}}$$

V is the potential difference. The magnetic field acting on the electron is a sum of the external magnetic field from outside sources and the field due to the coils. Thus the total field acting on the electron is:

$$B = B_C + B_e = \left(\frac{4}{5}\right)^{\frac{3}{2}} \frac{\mu_0 nI}{R} + \frac{\sqrt{2}R}{\mu_0 n} \left(\frac{5}{4}\right)^{\frac{3}{2}}$$

 I_0 is a constant proportional to the external magnetic field, μ_0 is defined in the appendix as **Constant 1**, R is the radius of the coils and n is the number of turns in each coil. Thus, this equation can be subbed into the equation describing the curvature of the electron path to obtain

$$\frac{1}{r} = k \frac{I - I_0}{\sqrt{V}} \sqrt{\frac{e}{m}}$$

Equation 1: Electron path model

From this one can determine the experimental value of the charge to mass ratio of an electron.

2 Procedure and Methodology

The experimental apparatus as pictured in **Figure 1** consisted of a lightbulb containing hydrogen gas and an electron gun. The electrons were accelerated through an anode connected to a power source which varies between 0-300V and a voltmeter which was set to 1000V. This was connected in a circuit with two Helmholtz coils whose current was provided by an 8V power supply in series with an ammeter which was set to 3A and a rheostat. To determine the radius of one Helmholz coil, the diameter was measured using a ruler and the number of turns in the coils was indicated on the apparatus.




Figure 1: Apparatus Schematic

Experiment 1: The 0-300V power supply and then the 8V power supply were turned on. The measurements of the diameter of the paths first began by fixing the 0-300V power supply voltage to a value of $149.820V \pm 0.009$. The current was varied using the rheostat for ten values incremented by 0.1A from 1-2A, each measured using the ammeter. For each current measurement, the diameter of the electron path was measured using a self-illuminated scale and plastic reflector and recorded. This procedure was repeated another time for values within a similar range of each increment.

Experiment 2: Similarly, the same procedure outlined above was repeated, however the current was fixed to a value of $1.511 \pm 0.003A$ using the rheostat. The voltage was then varied by increments of 10V from 200 - 300V and measured using the voltmeter. The diameter of the electron beam path was measured and the procedure was repeated for a second time for values within a similar range of each increment.

Note all the measurements were recorded in a Google Sheet and the voltmeter and ammeter were Keysight 34461A multimeters.

3 Results

Note, all referenced functions, equations and constants can be found in the appendix.

To determine the experimental value of the external magnetic field, the model as implemented by **Function 3** was defined. An array of coil field values was first determined using **Equation 2**, passing in the varied current data. The radius of the two coils was determined using a ruler to be $31.1 \pm 0.1cm$. Since the markings on the ruler were a bit faded, an uncertainty of 0.05cm for either end of the measured section was used. The number of turns was indicated on the apparatus and μ_0 (Vacuum Permeability) is a constant defined in the appendix (**Constant 1**). **Function 3** (implementation of **Equation 6**) was passed into scipy's curve_fit function and was fitted to the coil field and current data. The magnitude of the external magnetic field was determined to be $0.00011 \pm 0.00002T$. The uncertainty was determined using curve_fit's covariance values since the diagonals of the covariance matrix represent the variability of each parameter and the square root of this was taken to calculate the standard error. **Figure 2** demonstrates the linear relationship between the coil's field and the current where the y-intercept represents the value of the external field. **Figure 3** is the corresponding residuals plot for **Figure 2**.




Figure 2: Graph of fit used to determine the coil field




Figure 3: Residuals of coil field fit

Equation 5. Figures 4 and 5 demonstrate the relationship between the current and the radius of the electron beam and the voltage and the radius of the electron beam respectively.

0.0775 0.0750 0.0750




Figure 4: Experiment 1 data plot of electron path radius vs current through the Helmholtz coils. Includes best fit curve, error bars and the theoretical model where the literature charge mass ratio was used




Figure 5: Experiment 2 data plot of electron path radius vs potential difference across the electron gun anode. Includes best fit curve, error bars and the theoretical model where the literature charge mass ratio was used




Figure 6: Experiment 1 residual plot




Figure 7: Experiment 2 residual plot

Note the current and voltage uncertainties were determined from the Keysight 34461A multimeter manual and were implemented using **Function 4**. An uncertainty of 3*0.05cm was used as the uncertainty of the diameters of the electron beam paths due to measurements being taken on either side and the middle section of the self-illuminated scale. The self-illuminated scale and plastic reflector were used instead of a traditional ruler to eliminate problems of parallax. When a measurer uses the scale, the scale and reflector work to "project" the scale onto the electron beam path so that it seems like the scale is as close to the beam as possible. Thus viewing the beam from a different angle will not cause much discrepancy in the measured values.

Lastly the reduced chi-squared values determined for the varied current model and varied voltage model were 36.47 and 0.97 respectively using **Function 4**.

4 Analysis and Discussion

Throughout the experiment the electrons consistently followed the predicted circular path. Consider $\frac{1}{r} = \frac{B}{\sqrt{V}} \sqrt{\frac{e}{2m}}$, rearranging the equation for r, $r = \frac{\sqrt{V}}{B} \sqrt{\frac{2m}{e}}$. As $V \to 0$ and $B \to \infty$, $r \to 0$ and thus the electron's circular path becomes arbitrarily small. It follows that measuring the radius of the electrons' path in this case becomes more difficult as the trajectory becomes smaller, introducing the likelihood of greater error within the measurements. One way to reduce error, and this applies to trajectories of any radius, is to use a computer program which takes images of the circular electron beam and uses a scale to obtain precise measurements of the radius instead of using the self-illuminated scale and plastic reflector. This would also help to eliminate/reduce the effects of parallax on the obtained measured values.

The magnitude of the final value determined for the external magnetic field was approximately $0.00011 \pm 0.00002T$ which is close to the ideal value of 0T for a perfect experimental set up where no external field exists. Although an external magnetic field was present, other ferromagnetic materials and objects which generate magnetic fields such as cellphones did not noticeably affect the trajectory of the electrons when held close to the bulb. The Earth's magnetic field is approximately $25000 - 65000nT < 0.00011 \pm 0.00002T$, thus the experimental value indicates other sources contributed to the external magnetic field such as the building and other devices/instruments ("Geomagnetism Frequently Asked Questions").

The final experimental values determined for the electron charge to mass ratios from the varied current and voltage data were $2500000000000000 \pm 970000000000/kg$ and $1800000000000 \pm 900000000C/kg$ and the reduced chi-squared values were 36.47 and 0.97 respectively. In comparison, the theoretical value rounded to two significant figures is 18000000000C/kg. The experimental value determined from the data collected using a fixed value of current agrees with the theoretical value but although the order of magnitude is correct, the theoretical value does not fall within the range of the experimental value determined from the data collected using a fixed value of voltage. However, during the experiment the current reading on the ammeter significantly fluctuated which likely resulted in greater measurement error. Error may have also resulted from the actual measurement process as a better placement of the self-illuminated scale and plastic reflector was determined when data was being collected for a fixed value of current. In regards to the reduced chisquared values determined, the value 0.97 corresponding to data collected using a fixed current indicates that the corresponding model used is a good fit for the experimental data and examining the line of best fit in Figure 5, it is close to the theoretical line. In contrast, the reduced chi-squared value of 36.47 corresponding to data collected using a fixed voltage indicates that the model is ill-fitting for the experimental data and examining the curve of best fit in Figure 4, it is quite far from the theoretical curve. As stated earlier, error likely came from the experimental measurements. Figure 6 and 7 are the residual plots for Figures 4 and 5 respectively.

5 Conclusion

Although the experimental value determined from the data collected using a fixed value of voltage deviates away from the theoretical electron charge to mass ratio, this was likely due to measurement errors. In contrast, the experimental value determined from the data collected using a fixed value of current agrees with the theoretical value and thus verifies that the charge to mass ratio of an electron can be quantified using Newton's second law of motion.

6 References

National Oceanic and Atmospheric Administration. (n.d.). Geomagnetism Frequently Asked Questions.

National Centers for Environmental Information. Retrieved November 1, 2022, from

https://www.ngdc.noaa.gov/geomag/faqgeom.shtml#:~: text=The%20Earth's%20magnetic%20field%20intensity. magnetic%20north%20and%20true%20north.

7 Appendix

Constants

$$\mu_0 = 4\pi \cdot 10^{-7} WbA^{-1} m^{-1}$$

Constant 1: Vacuum Permeability

$$k = \frac{1}{\sqrt{2}} \left(\frac{4}{5}\right)^{\frac{3}{2}} \frac{\mu_0 n}{R}$$

Constant 2: k constant

Equations

$$\frac{1}{r} = k \frac{I - I_0}{\sqrt{V}} \sqrt{\frac{e}{m}}$$

Equation 1: Electron path model

$$u(x) = \pm \left| \frac{xp_v}{100} + \frac{rp_r}{100} \right|$$

Equation 3: Multimeter uncertainty calculation formula

$$u(x^2) = 2xu(x)$$

Equation 5: Squared value error propagation

$$B_c = \left(\frac{4}{5}\right)^{\frac{3}{2}} \frac{\mu_0 nI}{R}$$

Equation 2: Coil magnetic field model

$$\chi_R^2 = \frac{1}{N-n} \sum_{i=1}^{N} \left(\frac{y_i - y(x_i)}{u(y_i)} \right)$$

Equation 4: Reduced chi-Squared Metric

$$B_c = \frac{1}{r} \sqrt{\frac{2m}{e}} \sqrt{V} - B_e$$

Equation 6: Relation between radius, background field and coil field

Python Functions

```
def variedCurrent(I, a):
    return 1 / (a * k * ((I - externalFitVariables[1] / k) / np.sqrt(149.820)))
Function 1: Electron path model with constant voltage set to 149.820 Volts
def variedVoltage(v, a):
    return 1 / (a * k * ((1.510946 - externalFitVariables[1] / k) / np.sqrt(v)))
Function 2: Electron path model with constant current set to 1.510946 Amps
def externalField(r, b, B):
    return b * np.sqrt(149.820) / r - B
      Function 3: Coil magnetic field model (implements Equation 6)
def multimeterUncertainty(value, valuePercentage, range, rangePercentage):
    return (value * valuePercentage + range * rangePercentage) / 100
   Function 4: Function used to calculate uncertainty for the multimeter
                       (implements Equation 3)
def reducedChiSquared(x, y, yerr, model, modelParams):
    return (
        / (len(y) - len(modelParams))
        * np.sum(((y - model(x, *modelParams)) / yerr) ** 2)
    )
```

Function 5: Function used to calculate reduced chi-squared values (implements Equation 4)

Raw data tables

Current through	Measured elec-
Helmholtz coil	tron path radius
(A)	(m)
1.0078±0.002	0.071±0.0008
1.015223±0.002	0.072±0.0008
1.10219±0.003	0.068±0.0008
1.10813±0.003	0.064±0.0008
1.199013±0.003	0.061±0.0008
1.205404±0.003	0.057±0.0008
1.29518±0.003	0.053±0.0008
1.303277±0.003	0.056±0.0008
1.40726±0.003	0.052±0.0008
1.415226±0.003	0.050 ± 0.0008
1.508114±0.003	0.046±0.0008
1.513418±0.003	0.050 ± 0.0008
1.600502±0.003	0.043±0.0008
1.61092±0.003	0.044±0.0008
1.70728±0.004	0.039±0.0008
1.711507±0.004	0.039±0.0008
1.794713±0.004	0.038±0.0008
1.818418±0.004	0.036±0.0008
1.900371±0.004	0.034±0.0008
1.923341±0.004	0.034±0.0008
1.993017±0.004	0.032±0.0008
1.993969±0.004	0.031±0.0008

Table 1: Experiment 1 raw data: Varied current with constant anode voltage of 149.820V

Electron gun an- ode voltage (V)	Measured electron path radius (m)
199.334±0.01	0.066±0.0008
200.944±0.01	0.065±0.0008
210.022±0.01	0.067±0.0008
210.484±0.01	0.066±0.0008
220.791±0.01	0.070±0.0008
220.961±0.01	0.068±0.0008
229.437±0.01	0.071±0.0008
230.539±0.01	0.070±0.0008
240.056±0.01	0.071±0.0008
240.284±0.01	0.072±0.0008
250.261±0.01	0.074±0.0008
250.449±0.01	0.073±0.0008
259.95±0.01	0.076±0.0008
260.929±0.01	0.075±0.0008
269.85±0.01	0.077±0.0008
270.277±0.01	0.077±0.0008
279.812±0.01	0.078±0.0008
281.619±0.01	0.078±0.0008
289.422±0.01	0.079±0.0008
291.487±0.01	0.080±0.0008
300.142±0.01	0.081±0.0008
301.939±0.01	0.082±0.0008

Table 2: Experiment 2 raw data: Varied voltage with constant coil current of 1.510946A

Sample Calculations

This calculation is based on Equation 3

 $p_V = 0.002$, r = 1000, $p_r = 0.0006$ (from the Keysight manual)

Let V_i represent an arbitrary voltage measurement made with the 1000V setting.

$$u(V_i) = \pm \left| \frac{0.002V_i}{100} + \frac{1000 \times 0.0006}{100} \right| = \pm \left| 2 \times 10^{-5} V_i + 0.006 \right|$$

For example in table, the first measured voltage was: $V_1 = 199.334V$

$$u(V_1) = \pm |2 \times 10^{-5} \times 199.334V + 0.006| \approx \pm 0.00999$$

Calculation 1: Voltmeter Calculation sample for 1000V setting

This calculation is based on **Equation 3**:

 $p_V = 0.18$, r = 3, $p_r = 0.02$ (from the Keysight manual)

Let I_i represent an arbitrary current measurement made with the 3A setting.

$$u(I_i) = \pm \left| \frac{0.18I_i}{100} + \frac{3 \times 0.02}{100} \right| = \pm \left| 1.8 \times 10^{-3} I_i + 0.0006 \right|$$

For example in table one, the first measured current was: $I_1 = 1.0078A$.

$$u(I_1) = \pm |1.8 \times 10^{-3} \times 1.0078 + 0.0006| \approx \pm 0.0024$$

Calculation 2: Ammeter Calculation sample for 3A setting

In practice all calculations were done using the python implementation of **Equation 3**, **Function 4**.